 Medicom SAS
 Bld de la chanterie | 49124 Saint Barthélemy d'Anjou
 [BP 10059 – 49181 Saint Barthélemy d'Anjou Cedex] France
 Tel +33 0 241 963 434 | sales.support.pro.ue@medicom.com | www.medicom-eu.com
 Swiss Authorised Representative (CH-REP)
 FINAGEST SA, Rue Pierre-Fatio 12, 1204 Geneva - Switzerland
 Swiss Single Registration Number (CHRN): CHRN-AR-20002401

SafeSeal® Quattro


EN SafeSeal Quattro Self-Sealing Sterilization Pouches
 FR Poches de stérilisation auto-scellantes SafeSeal Quattro
 DE SafeSeal Quattro Selbstklebende Sterilisationsbeutel
 ES Bolsas de esterilización autosellantes SafeSeal Quattro
 IT SafeSeal Quattro Buste di sterilizzazione autosigillanti
 PL Samozamykające się torebki do sterylizacji SafeSeal Quattro
 NL SafeSeal Quattro zelfsluitende sterilisatiezakken
 PT Bolsas de esterilização auto-selantes SafeSeal Quattro
 UA Самоклеючі пакети для стерилізації SafeSeal Quattro

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EN INSTRUCTIONS FOR USE

INTENDED USE

Self-sealing sterilization pouch, are intended for use as packaging material for reusable instruments in sterilization by steam and ethylene oxide (ETO) in healthcare settings. The products are for single-use only. Products are designed to be used by trained healthcare professionals only.

LIMITATIONS FOR USE

Self-sealing sterilization pouches are not suitable for sterilization by irradiation, by hydrogen peroxide, formaldehyde, by hot, dry air or at temperature over 140 °C. Do not use pouch in dry heat sterilizers.

PACK INSTRUMENT IN THE POUCH - illustration 1

Select the correct pouch size for your sterilization items and ensure that the sterilization pouch is only filled to three quarters of its volume. Leave at least 3 cm of space around each instrument. Instrument should be placed in the pouch and not dropped so that the bottom seals remain intact. Medical device should be placed in the pouch in such a way that it can be easily removed to maintain sterility. Please note that sharp objects can damage the pouch and sterilization is therefore not guaranteed.

CLOSING THE POUCH - illustrations 2&3

All seals, including the closure, must be smooth and free from folds, bubbles or creases. The plastic film is darker where the seal has been applied. This makes it easier to check whether the seal is intact. No fasteners should be used to compress the packaging (e.g. rubber bands, paper clips, etc). Peel of protective strip. Fold over paper flap along fold and close the pouch by pressing the adhesive strap onto the transparent film. Press down from center outwards. The adhesive strap must be folded and glued precisely at the film to ensure a complete sealing.

AUTOClave

Slide the pouches vertically into autoclave or when filling horizontally, ensure that the pouches do not touch each other so that the entire surface of the pouch is exposed to steam or ethylene oxide. The pouches must not come into contact with the chamber walls. If a sterilization cycle needs to be repeated due to malfunction or a cycle is aborted before completion, all pouches must be replaced and repacked before being used in another sterilization cycle.

Follow the instructions of your device manufacturer for sterilization process.

AFTER STERILIZATION

The packaging and product must cool down before handling. Check each pouch to ensure that it has not been damaged during the procedure. Damaged pouches must be discarded and the devices reprocessed. Check the indicator has changed to the appropriate colour (steam/ETO). The bags must be clean and dry. If a pouch is wet, the devices inside should not be considered sterile. Place them in a new bag and process them again. Store items kept sterile in a dry closed place before use.

HOW TO VALIDATE THE STERILIZATION:

Chemical indicators are printed on each pouch (top and bottom corner). They change colour during the sterilization process as visual indication that the pouches have been sterilized -

illustrations 4

OPENING AND ASEPTIC PRESENTATION

The seals on the side with the process indicators should be unattached first. The package should then be opened by pulling the laminate away slowly and evenly. Open only 1/3 - 2/3 of the packaging length. When opening large and / or heavy packages they need to be supported by a table or a tray. Assistance may be needed to prevent any contamination of the packed instrument / instrument sets by accidentally touching the non-sterile outer surface of the packaging material.

SHELF LIFE OF PRODUCT

The product have a shelf of 3 years from date of manufacture. The expiry date is indicated on the packaging. Before initial use, always check that the product is within the stated shelf life. After steam or ETO sterilization, the pouch can be stored for a maximum 1 year while maintaining the sterility of the device until the sterile barrier is broken. This can be caused, for example, by tears, punctures, moisture or opening of the pouches.

INCIDENT NOTIFICATION

If any serious incident has occurred in relation to the use of device it should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

STORAGE AND HANDLING CONDITIONS

Sterilization pouches are single use disposable device and should not be reused. The device must be inspected before use: any defects, tears or damage also mean the pouch should be discarded and not used for sterilization. Do not use if package is damaged. Keep away from sunlight. Keep dry.

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