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EU declaration of conformity can be accessed is on: www.medicom-eu.com

Rev01 - 25/04/2024

MODEL OVERBOOTS: WL-CBE/AD-XX  
OVERSHOES: WL-CCE/AD-XX

CE 2895 Conformity to regulation EU - 2016/425  
Shirley EC Notified Body number 2895-Port Tunnel Business Park, Office 13 Unit 21, Dublin 17, ROI Ireland

02 ENGLISH	09 БЪЛГАРСКИ
02 FRANÇAIS	09 SLOVENŠČINA
03 ESPAÑOL	10 NORSK
04 DEUTSCH	10 LIETUVIŲ K.
04 ITALIANO	11 LATVIEŠU
05 NEDERLANDS	12 MAGYAR
05 PORTUGUÊS	12 SUOMI
06 POLSKI	13 EESTI
06 ROMÂNĂ	13 HRVATSKI
07 TÜRKÇE	14 SLOVENSKY
08 ÇEKSKY	14 ΕΛΛΗΝΙΚΑ
08 DANSK	15 СРПСКИ

ENGLISH  
1. Trademark  
2. Model identification  
3. CE marking : PPE category III  
4. European Protection Standards  
5. Refer to instructions for use  
6. Single use  
7. Fabric care symbols

**PHYSICAL PERFORMANCE OF WEEPRO ACCESSORIES®**

Fabric Physical Properties	Test Method	Class
Abrasion resistance	EN 530	1
Flex cracking	ISO 7854 B	4
Trapezoidal Tear MD / XD	ISO 9073-4	2
Tensile strength MD / XD	ISO 13934-1	1
Puncture resistance	EN 863	1
Seam strength	ISO 13935-2	3

**RESISTANCE TO PENETRATION BY LIQUIDS EN ISO 6530**

Fabric Chemical Properties	Penetration	Repellency
Sulphuric Acid 30 %	CLASS 3	CLASS 3
Sodium Hydroxide 10 %	CLASS 3	CLASS 3
Butan-1-ol	CLASS 3	CLASS 3
o-xylene	CLASS 3	CLASS 2

**FABRIC RESISTANCE TO PENETRATION OF INFECTIVE AGENTS EN14126**

Test Method	Class
Resistance to penetration by blood-borne pathogen ISO 16604	6/6
Resistance to penetration by contaminated liquids EN ISO 22610	6/6
Resistance to penetration by biologically contaminated aerosols EN ISO 22611	3/3
Resistance to penetration by contaminated solid particles ISO 22612	3/3

Typical Area of use  
Accessories are designed to protect against chemical risks (splashes – light sprays, Type PB [6]) and against infectious agents.

Limitations of use  
Accessories are intended to protect a partial body. The user will be the sole judge of the suitability of the protective accessories for the situation, the protection required, and the choice and combination of accessories with other equipment (gloves, respiratory protective equipment, etc.). Type PB [6] partial body protection has not been tested to the whole suit test according to EN 13034:2005+A1:2009. Accessories are for single use only and must be disposed after any job. Wear for long periods of time can cause heat stress. Heat stress and discomfort can be reduced or eliminated by using appropriate undergarments or suitable ventilation equipment. The manufacturer cannot accept responsibility for any improper use or disposal of garments produced by them. These accessories are flammable, keep away from fire. Do not use the accessories if any defects are noticed (e.g. holes, damaged seams, ...). Use the product within 5 years of the date of manufacturing.

Storage and disposal  
The accessories must be stored in their original packaging in a dry place, away from sunlight. They must be disposed of in compliance with current regulations.

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Resistances to penetration by blood-borne pathogen ISO 16604  
Resistance to penetration by contaminated liquids EN ISO 22610  
Resistance to penetration by biologically contaminated aerosols EN ISO 22611  
Resistance to penetration by contaminated solid particles ISO 22612

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