ENGLISH

1. Trademark	5. Size and body measurements
2. Reference	6. Refer to instructions for use
3. CE marking : PPE Category III	7. Fabric care symbols
4. European Protection Standards	8. Single use

PHYSICAL PERFORMANCE OF WEEPRO MAX GREEN®

Fabric Physical Properties	Test Method	Class
Abrasion resistance	EN 530	1
Flex cracking resistance	ISO 7854 B	4
Trapezoidal Tear MD / XD	ISO 9073-4	2
Tensile strength	ISO 13934-1	1
Puncture resistance	EN 863	1
Seam strength	ISO 13935-2	3

RESISTANCE TO CHEMICAL LIQUID EN ISO 6530

Fabric Chemical Properties	Penetration	Repellency
Sulphuric Acid 30%	Class 3	Class 3
Sodium Hydroxide 10%	Class 3	Class 3

FABRIC AND TAPE SEAMS RESISTANCE TO PERMEATION BY LIQUIDS EN ISO 6529

Fabric Chemical Properties	Breakthrough time (min)	Class
Sulphuric Acid 30%	10-30 min	1/6

FABRIC RESISTANCE TO PENETRATION OF INFECTIVE AGENTS EN14126

Test Method	Class
Resistance to penetration by blood-borne pathogens ISO 16604	6/6
Resistance to penetration by contaminated liquids EN ISO 22610	6/6
Resistance to penetration by contaminated liquid aerosols EN ISO 22611	3/3
Resistance to penetration by contaminated solid particles ISO 22612	3/3

WHOLE SUIT TEST PERFORMANCE

High level spray test EN 14605	Pass
Inward leakage test Method as defined by EN ISO 13982-1:2004 + A1:2010	Pass
Light splashes Method as defined by EN 13034 - 2005 + A1: 2009	Pass
Nominal protection factor according to EN 1073-2 - 2002	Class 2
Surface resistivity EN1149-5 - 2018	Pass

PROTECTION AGAINST PESTICIDES

Tests	Test Method	Requirements level C2	Results
Material resistance to penetration	ISO 22608	≤ 5 %	2.3 %
Material repellency	ISO 22608	≥ 80 %	93.9 %
Seam penetration resistance	ISO 22608	≤ 40%	0.0 %
Practical performance test	EN ISO 27065	Pass	Pass

Typical Area of use

Green coverall with elasticated hood, wrists, ankles and back, adhesive flap over the zipper, thumb loop. Coverall is designed to protect against chemical risks (splashes – sprays – airborne solid particulates) and against infectious agents and against contamination by radioactive particles and electrostatics (dissipative properties). The coverall protects against the use of insecticides, herbicides, fungicides, and other substances applied in liquid form (level C2).

How to wear protective clothing

Remove the coverall from its packaging, open the center zipper and put on the garment, taking care not to damage it. Close the zipper completely and cover it with the adhesive flap, without creasing. If there is a risk of airborne solid particles, it is advisable to tape additionally protective gloves with ends of the sleeves.

Limitations of use

This protective clothing is not suitable for use with concentrated pesticide formulations and against fumigants. The coverall do not protect against biocidal products used for agricultural and non-agricultural settings. Do not use the coverall if any defects are noticed (e.g., holes, damaged seams and fastenings, ...). The coverall meets the requirements: Ljm 82/90 S30% - Ls 8/10 ≤ 15%. Type 6 chemical protective suits have been tested to the whole suit test. Wear for long periods of time can cause heat stress. Heat stress and discomfort can be reduced or eliminated by using appropriate undergarments or suitable ventilation equipment. Coverall is for single use only and must be disposed after any job. The user is the sole judge of the adequacy between his risk analysis and the type of protection required, including the additional equipment necessary for a complete protection of the body. If necessary, use additional PPE with the same level of protection (such as as gloves, breathing apparatus, etc.) in order to provide full body protection. The manufacturer cannot accept responsibility for any improper use or disposal of garments produced by them. This garment is flammable, keep away from fire. The person wearing the electrostatic dissipative protective clothing shall be properly earthed. The resistance between the person skin and earth shall be less than 10⁶ Ω e.g. by wearing adequate footwear on dissipative or conductive floors. Electrostatic dissipative protective clothing shall not be open or removed whilst in presence of flammable or explosive atmospheres or while handling flammable or explosive substances. Electrostatic dissipative protective clothing is intended to be worn in Zones 1, 2, 20, 21 and 22 (see EN 60079-10-1 and EN 60079-10-2) in which the minimum ignition energy of any explosive atmosphere is not less than 0.016 mJ. Electrostatic dissipative protective clothing shall not be used in oxygen enriched atmospheres, or in Zone O (see EN 60079-10-1) without prior approval of the responsible safety engineer. The electrostatic dissipative performance of the electrostatic dissipative protective clothing can be affected by wear and tear, laundering and possible contamination electrostatic dissipative protective clothing shall be worn in such a way that it permanently covers all non-compliant materials during normal use. Use the product within 5 years of the date of manufacturing.

Storage and disposal

Coveralls should be transported and stored in the original packaging in a dry place, away from sunlight. They must be disposed of in compliance with current regulations.

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