

**1121**

## LATEX GLOVE SAFETOUCH VITALS POWDERED

Class I Medical Device and Personal Protective Equipment Category III

### FEATURES & BENEFITS

- Great elasticity.
- Good dexterity.
- Watertight and airtight.
- Contraindication: Latex is a natural material that contains protein which might cause allergies.

### Intended Use:

- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from tested chemicals and microorganisms (bacteria, moulds and viruses).




### PRODUCT DESCRIPTION

Product name :	<b>SafeTouch Vitals</b>
Product type :	Single use, non-sterile
Material :	Latex
Colour :	Natural
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured fingertips
Powdered :	Yes
Length :	240 mm
Unit Weight (M):	5g ±0,2g
Palm thickness :	min 0,08 mm
Finger thickness :	min 0,08 mm
Country of origin :	Thailand

Item	Size	
1121-A	T 6	XS
1121-B	T 7	S
1121-C	T 8	M
1121-D	T 9	L
1121-E	T 10	XL

Marking on packaging:



### INDUSTRIES

Health / Medical / Dental / Industry

## CERTIFICATION AND STANDARDS

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.

Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.

## PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

For single use only, do not reuse.

Please refer to your national legislation about disposal of this product once it has been used.

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at room temperature.

Shelf life of the product: 5 years.

## LOGISTIC INFORMATION








### Specifications outer case



### Specifications inner case

### Colour code dispenser and carton

Item	Size mm	Gross weight kg	QTY/ pallet	Size mm	QTY	Colour
1121-A	225 x 310 x 210	4,5	60	200 x 110 x 60	10 x 100	
1121-B	225 x 310 x 210	5,1	60	200 x 110 x 60	10 x 100	
1121-C	225 x 310 x 210	5,6	60	200 x 110 x 60	10 x 100	
1121-D	225 x 310 x 210	6,1	60	200 x 110 x 60	10 x 100	
1121-E	225 x 310 x 210	6,6	60	200 x 110 x 60	10 x 100	



### EAN Codes

Item	Inner case	Outer case
1121-A	3 662 036 021 489	3 662 036 021 533
1121-B	3 662 036 021 496	3 662 036 021 540
1121-C	3 662 036 021 502	3 662 036 021 557
1121-D	3 662 036 021 519	3 662 036 021 564
1121-E	3 662 036 021 526	3 662 036 021 571