

1122

LATEX GLOVE SAFETOUCH VITALS POWDER-FREE

Class I Medical Device and Personal Protective Equipment Category III

FEATURES & BENEFITS

- Great elasticity.
- Good dexterity.
- Watertight and airtight.
- Contraindication: Latex is a natural material that contains protein which might cause allergies.



Intended Use:


- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from tested chemicals and microorganisms (bacteria, moulds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.
- For food handling / preparation (HACCP).

PRODUCT DESCRIPTION

Product name :	SafeTouch Vitals
Product type :	Single use, non-sterile
Material :	Latex
Colour :	Natural
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Unit Weight (M):	5,5g ±0,2g
Palm thickness :	min 0,08 mm
Finger thickness :	min 0,08 mm
Country of origin :	Thailand

Item	Size	
1122-A	T 6	XS
1122-B	T 7	S
1122-C	T 8	M
1122-D	T 9	L
1122-E	T 10	XL

Marking on packaging:



INDUSTRIES

Health / Medical / Dental / Industry/ Food

CERTIFICATION AND STANDARDS

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.

Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.

In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August, 2020.

This glove is suitable for food contact (restriction: do not use with acidic food).

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

For single use only, do not reuse.

Please refer to your national legislation about disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at room temperature.

Shelf life of the product: 5 years.

LOGISTIC INFORMATION








Specifications outer case



Specifications inner case

Colour code dispenser and carton

Item	Size mm	Gross weight kg	QTY/ pallet	Size mm	QTY	Colour
1122-A	249 x 340 x 250	4,5	60	240 x 122 x 65	10 x 100	
1122-B	249 x 340 x 250	5,1	60	240 x 122 x 65	10 x 100	
1122-C	249 x 340 x 250	5,6	60	240 x 122 x 65	10 x 100	
1122-D	249 x 340 x 250	6,1	60	240 x 122 x 65	10 x 100	
1122-E	249 x 340 x 250	6,6	60	240 x 122 x 65	10 x 100	



EAN Codes

Item	Inner case	Outer case
1122-A	3 662 036 021 588	3 662 036 021 632
1122-B	3 662 036 021 595	3 662 036 021 649
1122-C	3 662 036 021 601	3 662 036 021 656
1122-D	3 662 036 021 618	3 662 036 021 663
1122-E	3 662 036 021 625	3 662 036 021 670