

1130

LATEX GLOVE SAFETOUCH CONNECT SOFT POWDER-FREE

Class I Medical Device and Personal Protective Equipment Category III

FEATURES & BENEFITS

- Great elasticity and good dexterity.
- Faster gloving thanks to its internal polymer coating.
- Watertight and airtight.
- Contraindication: Latex is a natural material that contains protein which might cause allergies.



Intended Use:

- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from tested chemicals and microorganisms (bacteria, moulds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.

PRODUCT DESCRIPTION

Product name :	SafeTouch Connect Soft
Product type :	Single use, non-sterile
Material :	Latex
Colour :	Natural
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Unit Weight (M):	6,2g ± 0,2g
Palm thickness :	min 0,08 mm
Finger thickness :	min 0,08 mm
Country of origin :	Thailand

Item	Size	
1130-A	T 6	XS
1130-B	T 7	S
1130-C	T 8	M
1130-D	T 9	L
1130-E	T 10	XL

Marking on packaging:



INDUSTRIES

Health / Medical / Dental / Industry

CERTIFICATION AND STANDARDS

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, UK Medical Device Regulation 2002 and the Swiss Ordonnance (ODim) from July 1st 2020, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425 and Regulation 2016/425 as amended to apply in GB, on personal protective equipments.

Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

For single use only, do not reuse.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at room temperature.

Shelf life of the product: 5 years.

LOGISTIC INFORMATION








Specifications outer case



Specifications inner case

Colour code dispenser and carton

Item	Size mm	Gross weight kg	QTY/ pallet	Size mm	QTY	Colour
1130-A	249 x 340 x 250	6,3	63	240 x 122 x 65	10 x 100	
1130-B	249 x 340 x 250	6,8	63	240 x 122 x 65	10 x 100	
1130-C	249 x 340 x 250	7,5	63	240 x 122 x 65	10 x 100	
1130-D	249 x 340 x 250	7,9	63	240 x 122 x 65	10 x 100	
1130-E	249 x 340 x 250	8,2	63	240 x 122 x 65	10 x 100	



EAN Codes

Item	Inner case	Outer case
1130-A	3 662 036 021 380	3 662 036 021 434
1130-B	3 662 036 021 397	3 662 036 021 441
1130-C	3 662 036 021 403	3 662 036 021 458
1130-D	3 662 036 021 410	3 662 036 021 465
1130-E	3 662 036 021 427	3 662 036 021 472