

1176

NITRILE GLOVE SAFETOUCH REJUVENATE POWDER-FREE

Enriched in Lanolin and Vitamin E

Class I Medical Device and Category III Personal Protective Equipment

FEATURES & BENEFITS

- Great elasticity.
- Good dexterity.
- Its micro-textured appearance gives it correct touch and grip of objects.
- Enriched in Lanolin and Vitamin E, this glove gives the feeling of a better hydrated skin, softer and smoother.



Recommendation:

- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from microorganisms (bacteria, moulds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.

PRODUCT DESCRIPTION

Product name :	SafeTouch Rejuvenate
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	White
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Fingertips
Powdered :	No
Length :	240 mm
Palm thickness :	0,11±0,03 mm
Finger thickness :	0,13±0,03 mm
Integrated lotion:	Lanolin and Vitamin E
Unit weight (M):	4g (±0,2g)
0,2g)Country of origin :	Malaysia

Item	Size	
1176-A	T 6	XS
1176-B	T 7	S
1176-C	T 8	M
1176-D	T 9	L
1176-E	T10	XL

Marking on packaging:



INDUSTRIES

Health / Medical/ Dental

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.

Comply with the EN ISO 21420, EN ISO 374-2/5, ISO 16604 for protective gloves against microorganisms including viruses.

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

Single use, do not reuse.

Please refer to your national legislation about disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.

Shelf life of the product: 5 years.

LOGISTIC INFORMATION



Specifications outer case

Item	Size mm	Gross weight kg	QTY/ Pallet
1176-A	320 x 265 x 250	4,4	63
1176-B	320 x 265 x 250	4,7	63
1176-C	320 x 265 x 250	5,1	63
1176-D	320 x 265 x 250	5,3	63
1176-E	320 x 265 x 250	5,5	63



Specifications inner case

Size mm	QTY
235 x 130 x 60	10 x 100
235 x 130 x 60	10 x 100
235 x 130 x 60	10 x 100
235 x 130 x 60	10 x 100
235 x 130 x 60	10 x 100

Colour code inner and outer

Colour



EAN Codes

Item	Inner case	Outer case
1176-A	3 700 381 514 802	3 700 381 514 413
1176-B	3 700 381 514 819	3 700 381 514 420
1176-C	3 700 381 514 826	3 700 381 514 437
1176-D	3 700 381 514 833	3 700 381 514 444
1176-E	3 700 381 521 565	3 700 381 521 558