

# 1191N

## NITRILE GLOVE SAFETOUCH VITALS POWDER-FREE

Class I Medical Device and Category III Personal Protective Equipment

### FEATURES & BENEFITS

- Latex-like feel.
- Agreeable dexterity.
- Excellent grip in both wet and dry conditions.
- Watertight and airtight.
- Convenient packaging that holds **200 gloves** per box reducing storage space, transportation costs and packaging waste.



### RECOMMENDATION:

- Users allergic to latex.
- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the healthcare professional, the patient and their environment.
- In an industrial environment, this glove is designed to protect the wearer from tested chemicals and microorganisms (bacteria, molds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.
- For food handling / preparations (HACCP).

### PRODUCT DESCRIPTION

Product name:	<b>SafeTouch® Vitals</b>
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	Blue
Shape :	Ambidextrous
Cuff:	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Palm thickness :	0,07 mm (±0,02mm)
Finger thickness :	0,09 mm (±0,02mm)
Unit weight (M):	3,2 g

Item	Size	
1191N-A	T 6	XS
1191N-B	T 7	S
1191N-C	T 8	M
1191N-D	T 9	L
1191N-E	T 10	XL

Marking on packaging:



### INDUSTRIES

Food / Health / Medical / Dental/ Administration / Hygiene/ Industry

## CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.

Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.

In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August 2020. This glove is suitable for food contact, with all types of foodstuff.

## PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation about disposal of this product once it has been used.

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored in a cool environment or at room temperature.

Shelf life of the product: 5 years.

## LOGISTIC INFORMATION








### Specifications outer case



### Specifications inner case

### Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/ pallet	Size mm	QTY	Colour
1191N-A	415 x 261 x 260	6,2	49	243 x 125 x 81	10 x 200	
1191N-B	415 x 261 x 260	6,5	49	243 x 125 x 81	10 x 200	
1191N-C	415 x 261 x 260	7,1	49	243 x 125 x 81	10 x 200	
1191N-D	415 x 261 x 260	8,4	49	243 x 125 x 81	10 x 200	
1191N-E	415 x 261 x 260	10,3	49	243 x 125 x 81	10 x 200	



### EAN Codes

Item	Inner case	Outer case
1191N-A	3 700 381 519 890	3 700 381 519 883
1191N-B	3 700 381 519 913	3 700 381 519 906
1191N-C	3 700 381 519 937	3 700 381 519 920
1191N-D	3 700 381 519 951	3 700 381 519 944
1191N-E	3 700 381 519 975	3 700 381 519 968