

1193N

NITRILE GLOVE SAFETOUCH VITALS POWDER-FREE

Class I Medical Device and Category III Personal Protective Equipment

FEATURES & BENEFITS

- Latex-like feel.
- Agreeable dexterity.
- Excellent grip in both wet and dry conditions.
- Watertight and airtight.

RECOMMENDATION:

- Users allergic to latex.
- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the healthcare professional, the patient and their environment.
- In an industrial environment, this glove is designed to protect the wearer from tested chemicals and microorganisms (bacteria, molds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.
- For food handling / preparations (HACCP).



PRODUCT DESCRIPTION

Product name:	SafeTouch® Vitals
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	Blue
Shape :	Ambidextrous
Cuff:	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Palm thickness :	0,07 mm (±0,02mm)
Finger thickness :	0,09 mm (±0,02mm)
Unit weight (M):	3,2 g

Item	Size	
1193N-A	T 6	XS
1193N-B	T 7	S
1193N-C	T 8	M
1193N-D	T 9	L
1193N-E	T 10	XL

Marking on packaging:



INDUSTRIES

Food / Health / Medical / Dental/ Administration / Hygiene/ Industry

1193N

NITRILE GLOVE SAFETOUCH VITALS POWDER-FREE

Class I Medical Device and Category III Personal Protective Equipment

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.
In accordance with requirements of Reg (EU) 2017/745, on medical devices.
Comply with the standards EN 455-1/2/3/4 for single use medical gloves.
In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.
Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.
In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August 2020. This glove is suitable for food contact, with all types of foodstuff.

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.
The reuse or extended use may lead to infection or cross-contamination.
Please refer to your national legislation about disposal of this product once it has been used

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored in a cool environment or at room temperature.
Shelf life of the product: 5 years.

LOGISTIC INFORMATION



Specifications
outer case



Specifications
inner case

Colour code
inner and outer

Item	Size mm	Gross weight kg	QTY/ pallet	Size mm	QTY	Colour
1193N-A	262 x 252 x 270	3,2	60	245 x 125 x 48	10 x 100	
1193N-B	262 x 252 x 270	3,4	60	245 x 125 x 48	10 x 100	
1193N-C	262 x 252 x 270	3,7	60	245 x 125 x 48	10 x 100	
1193N-D	262 x 252 x 270	4,5	60	245 x 125 x 48	10 x 100	
1193N-E	262 x 252 x 270	5,1	60	245 x 125 x 48	10 x 100	



EAN Codes

Item	Inner case	Outer case
1193N-A	3 700 381 519 999	3 700 381 519 982
1193N-B	3 700 381 520 018	3 700 381 520 001
1193N-C	3 700 381 520 032	3 700 381 520 025
1193N-D	3 700 381 520 056	3 700 381 520 049
1193N-E	3 700 381 520 070	3 700 381 520 063