

40.466

SPUNLACE BOXER SHORT

Cat 1 Personal Protective Equipment.

FEATURES & BENEFITS


- The material made from viscose ensures comfort for the wearer. It is soft, does not scratch, and is not irritating. It also absorbs sweat.
- White and opaque, it ensures the wearer's privacy.
- Loose fit for ease of movement.
- Elastic waistband.

PRODUCT DETAILS

Product name :	Spunlace boxer short
Product type :	Single use, non sterile
Material :	Viscose and polyester
Confection :	Sewn and elasticated at the waist
Grammage :	60 g/m ²
Unit weight :	35 g (+/- 3%)
Color:	White
Qty per case :	3 x 100 units
Origin :	China



Dimensions (+/- 1cm) :



CERTIFICATIONS, STANDARDS & TESTS

Complies with the requirements of European Regulation 2016/425 governing Personal Protective Equipment.

MECHANICAL RESISTANCE			
Test	Test Method	Test Report N°	Result
Tensile strength	ISO 9073-3	2401-032-5	Machine direction : 76,0 N Cross direction : 44,0 N
Elongation at maximum force	ISO 9073-3	2401-032-5	Machine direction : 33,0 % Cross direction : 73,4 %

40.466

SPUNLACE BOXER SHORT

Cat 1 Personal Protective Equipment.

STORAGE CONDITIONS

Normal storage and conservation conditions: should not be exposed to moisture and sunlight, should be stored at room temperature.

Product shelf life: 5 years

LOGISTICS

Item	Qty per case	Case size (cm)	Gorss weight (kg)	Case qty per pallet	Pallet height	EAN
40.466	3 x 100 pcs	50 x 38 x 48	11,7	12 cases	1,64 m	3 700 381 522 838

