

DISPOSABLE FACE DRAPE

Medical Device, Class I

TECHNICAL DATA SHEET

English

FEATURES & BENEFITS

- Consist of a net in the upper area as well as a lower part made of spunbonded polypropylene with a pre-punched hole
- Will be placed on the patient's face as shown in figure 1 to maintain visibility and stop patient from being covered in splashes
- Pre-punched hole allows good access to the patient's mouth
- Single-use



PRODUCT DESCRIPTION

Product name: Medicom SafeDrape™ disposable face drapes

Product type: Single use, non sterile

Material: Spunbond polypropylene / Nylon thread texturized polyester

Size: 380 mm (±10 mm) x 340 mm (±20 mm)

Packaging: Box 50 units

Origin: China

Available in 2 colours:

	White		Blue
	8-1110W		8-1110B

For professional use only.

INDUSTRIES

Medical / Dental

CERTIFICATION & STANDARD

Meet the requirements of Regulation (EU) 2017/745 on medical devices.

Manufacturing site ISO 13485 certified.



DISPOSABLE FACE DRAPE

Medical Device, Class I

TECHNICAL DATA SHEET

English

STORAGE CONDITIONS

Normal conditions of conservation and storage:
not to be exposed to moisture and sunlight.
The shelf life for disposable face drape is 2 years.

LOGISTICAL INFORMATION

Outer case specifications



REF	Size (mm)	EAN Code
8-1110B	415 x 200 x 410	3662036023353
8-1110W	415 x 200 x 410	3662036023377

Inner case (box) specifications



REF	Size (mm)	QTY/ Outer case	EAN Code
8-1110B	395 x 195 x 35	10 boxes x 50 pcs.	3662036023360
8-1110W	395 x 195 x 35	10 boxes x 50 pcs.	3662036023384