

## REUSABLE FACE DRAPE

Medical Device, Class I

## TECHNICAL DATA SHEET

English

### FEATURES & BENEFITS

- Consist of a net in the upper area as well as a lower part of cotton fiber towel with a hole
- Will be placed on the patient's face as shown in figure 1 to maintain visibility and stop patient from being covered in splashes
- Hole allows good access to the patient's mouth
- Washable



### PRODUCT DESCRIPTION

Product name: Medicom SafeDrape™ reusable face drapes

Product type: Washable

Material: Cotton fiber towel / Nylon thread texturized polyester

Size: 380 mm (±10 mm) x 340 mm (±10 mm)

Packaging: Box 12 units

Origin: China

Available in 5 versions:

	White 8-4020W		Blue 8-4020B		Yellow 8-4020G		Mint 8-4020M		Assorted* 8-4020S
-------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------	-----------------	-------------------------------------------------------------------------------------	-------------------	-------------------------------------------------------------------------------------	-----------------	-------------------------------------------------------------------------------------	----------------------

\* Assorted pack contains 4 each of blue, yellow and mint face drapes

For professional use only.


### INDUSTRIES

Medical / Dental

### CERTIFICATION & STANDARD

Meet the requirements of Regulation (EU) 2017/745 on medical devices.

Manufacturing site ISO 13485 certified.



## REUSABLE FACE DRAPE

Medical Device, Class I

## TECHNICAL DATA SHEET

English

### STORAGE CONDITIONS

Normal conditions of conservation and storage:  
not to be exposed to moisture and sunlight.

The shelf life for reusable face drape is 2 years.

### LOGISTICAL INFORMATION

#### Outer case specifications



REF	Size (mm)	EAN Code
8-4020W	415 x 250 x 410	3662036023476
8-4020B	415 x 250 x 410	3662036023391
8-4020G	415 x 250 x 410	3662036023414
8-4020M	415 x 250 x 410	3662036023438
8-4020S	415 x 250 x 410	3662036023452

#### Inner case (box) specifications



REF	Size (mm)	QTY/ Outer case	EAN Code
8-4020W	395 x195 x 45	10 boxes x 12 pcs.	3662036023483
8-4020B	395 x195 x 45	10 boxes x 12 pcs.	3662036023407
8-4020G	395 x195 x 45	10 boxes x 12 pcs.	3662036023421
8-4020M	395 x195 x 45	10 boxes x 12 pcs.	3662036023445
8-4020S	395 x195 x 45	10 boxes x 12 pcs.	3662036023469