



8280-N
to
8294-N

MEDICOM® DRY-BACK DENTAL BIBS 3-PLY

FEATURES & BENEFITS

- The Medicom® Dry-back 3-ply dental bibs range is made of 2 layers of embossed cellulose and 1 layer of waterproof polyethylene material.
- This allows maximum watertightness with horizontal embossed pattern and superior absorption thanks to additional layer of cellulose.

Recommended use :

This bib, placed polyethylene face on the user, permits to establish a watertight barrier against splashing liquid.



Available in 13 colours:

8280-N  Lavender	8282-N  Blue	8283-N  White	8284-N  Yellow	8285-N  Green	8286-N  Dusty Rose	
8288-N  Peach	8289-N  Aqua	8290-N  Dark blue	8291-N  Orange	8292-N  Black	8293-N  Lime	8294-N  Turquoise

PRODUCT DESCRIPTION

Product name:	Dry-Back Dental Bibs 3-Ply
Product type:	Single use, non-sterile
Material:	Embossed cellulose, Polyethylene
Colour :	Lavender, Blue, White, Yellow, Green, Dusty Rose, Peach, Aqua, Dark blue, Orange, Black, Lime, Turquoise
Grammage :	33 g/m ² (± 10%)
Folding:	M
Size:	33 x 45,5 cm (± 0,5 cm)
Unit weight:	6 g (± 10%)
Packaging:	Bulk of 500 units
Origin:	France



8280-N
to
8294-N

MEDICOM® DRY-BACK DENTAL BIBS 3-PLY

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C.

LOGISTICAL INFORMATION



Outer case specifications

Item	Size mm	Gross weight kg	QTY/ carton	QTY/ pallet
8280-N	300 x 240 x 340	3,2	Bulk of 500 pieces	78 cartons (6 layers of 13 cartons)
8282-N				
8283-N				
8284-N				
8285-N				
8286-N				
8288-N				
8289-N				
8290-N				
8291-N				
8292-N				
8293-N				
8294-N				



EAN Codes

Item	Outer case
8280-N	8 718 007 006 230
8282-N	8 718 007 006 377
8283-N	8 718 007 006 490
8284-N	8 718 007 006 605
8285-N	8 718 007 006 711
8286-N	8 718 007 006 810
8288-N	8 718 007 006 933
8289-N	8 718 007 007 039
8290-N	8 718 007 007 138
8291-N	8 718 007 007 183
8292-N	8 718 007 012 422
8293-N	8 718 007 013 719
8294-N	8 718 007 906 141



KOLMI HOPEN GROUP SAS

Bld de la chanterie | 49124 Saint Barthélemy d'Anjou – France

Bld de la chanterie | BP 10059 – 49181 Saint Barthélemy d'Anjou Cedex – France

T : + 33 (0) 241 963 434 | F : + 33 (0) 241 963 453 | sales@kolmi-hopen.com | www.medicom-eu.com