

**88000, 88005, 88010,
88015, 88025, 88030,
88035, 88040**


MEDICOM® SAFESEAL® QUATTRO STERILIZATION POUCHES

Medical Device, Class I according to MDR 2017/745

FEATURES & BENEFITS

- For use as sterile barrier packaging for terminally sterilized medical devices in sterilization by steam, ethylene oxide gas in healthcare establishments
- Non-tearing, thick multi-layer film made of a combination of polypropylene and polyethylene
- Transparent, blue tinted film facilitates both easy instrument viewing and identification of puncture or tears in the film
- Clean, fibre-free separation of paper and film ensures safe presentation of sterile products
- Pouches have a wide adhesive sealing tape which easily and securely seals the pouch providing a reliable seal until pouch is opened
- Permit to keep the sterilization conditions of the contained devices after sterilization: up to 1 year

Process Indicator Colour Change Guide:



PRODUCT DESCRIPTION

Product name: SafeSeal® Quattro Sterilization Pouches
 Product type: Single use, non-sterile
 Material: Heavy weight medical paper, polypropylene and polyethylene film
 Medical grade paper: 60 g/m²
 Film: 50 µm

Sterilization Cycle Parameter	Prevacuum Steam		Ethylene Oxide Gas
Time (min)	30	4	100- 120
Temperature (°C)	121	134	50
Relative humidity (%)	N/A	N/A	60 - 85
Sterilant Concentration	N/A	N/A	600
Drying-Aeration time	15 - 30 min	20 - 30 min	12 hrs

N/A = Not applicable

For color change guide, facts on the Stated Value (SV) and all other information, consult www.medicom.com

Origin: Taiwan

INDUSTRIES

Medical / Dental



REF	Size mm
88000	70 x 229
88005	89 x 133
88010	89 x 229
88015	57 x 102
88025	133 x 254
88030	191 x 330
88035	254 x 356
88040	305 x 432



**88000, 88005, 88010,
88015, 88025, 88030,
88035, 88040**

**MEDICOM® SAFESEAL® QUATTRO
STERILIZATION POUCHES**
Medical Device, Class I

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored at a temperature between 5 °C and 40 °C. Shelf life: 3 years.

LOGISTICAL INFORMATION



Outer case specifications

REF	Size mm	Gross weight kg	QTY/ pallet
88000	268 x 160 x 202	2,9	180
88005	483 x 173 x 74	2,1	200
88010	298 x 267 x 138	3,7	120
88015	328 x 141 x 138	2,2	285
88025	298 x 298 x 352	10,0	48
88030	328 x 210 x 370	8,5	60
88035	394 x 330 x 267	12,3	42
88040	476 x 337 x 343	17,7	25



Inner case specifications

Size mm	QTY/ Outer case
260 x 64 x 76	6 packages x 200 pcs.
165 x 64 x 95	5 packages x 200 pcs.
260 x 64 x 95	6 packages x 200 pcs.
133 x 64 x 64	10 packages x 200 pcs.
286 x 64 x 140	10 packages x 200 pcs.
362 x 64 x 197	5 packages x 200 pcs.
387 x 64 x 260	5 packages x 200 pcs.
464 x 64 x 311	5 packages x 200 pcs.



EAN Codes

REF	Inner case	Outer case
88000	8718007911978	8718007011975
88005	8718007911985	8718007011982
88010	8718007911992	8718007011999
88015	8718007912005	8718007012002
88025	8718007912012	8718007012019
88030	8718007912029	8718007012026
88035	8718007912036	8718007012033
88040	8718007912043	8718007012040



Medicom Healthcare B.V.
Parallelweg 80a | 3931 MT Woudenberg | The Netherlands
cs@medicom-eu.com | +31 332 866 460