

**9850, 9860, 9870, 9880,
9890, 98100, 98110**

**MEDICOM® SAFESEAL® QUATTRO
STERILIZATION REELS**


FEATURES & BENEFITS

- Guaranteed seal during sterilization
- Extremely easy to open without releasing fibres or tearing
- Offers maximum microbial protection
- Manufactured in full compliance with MDR and the EN 868-5 standard

PRODUCT DESCRIPTION

Product name: SafeSeal® Quattro Sterilization Reels
 Product type: Single use, non-sterile
 Material: Heavy weight medical paper and a double layer of light blue polyethylene/polypropylene film
 Grammage: 60 g/m² (Paper)

Item	Size mm
9850	50 mm x 200 m
9860	75 mm x 200 m
9870	100 mm x 200 m
9880	150 mm x 200 m
9890	200 mm x 200 m
98100	250 mm x 200 m
98110	300 mm x 200 m



Origin: Made in EU

INDUSTRIES

Medical / Dental

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 8 °C and 35 °C.

9850, 9860, 9870, 9880,
9890, 98100, 98110

MEDICOM® SAFESEAL® QUATTRO STERILIZATION REELS

LOGISTICAL INFORMATION



Outer case specifications

Content

Item	Size mm	Gross weight kg	QTY/ pallet	QTY/ Outer case
9850	225 x 225 x 420	8,90	60 cartons	8 Sterilization rolls
9860	225 x 225 x 420	8,40	60 cartons	5 Sterilization rolls
9870	225 x 225 x 420	9,10	60 cartons	4 Sterilization rolls
9880	225 x 225 x 420	6,90	60 cartons	2 Sterilization rolls
9890	225 x 225 x 420	9,10	60 cartons	2 Sterilization rolls
98100	225 x 225 x 420	5,80	60 cartons	1 Sterilization roll
98110	225 x 225 x 420	6,75	60 cartons	1 Sterilization roll



EAN Codes

Item	Outer case	Label
9850	8720249901952	8720249902027
9860	8720249901969	8720249902034
9870	8720249901976	8720249902041
9880	8720249901983	8720249902058
9890	8720249901990	8720249902065
98100	8720249902003	8720249902072
98110	8720249902010	8720249902089



Medicom Healthcare B.V.
Parallelweg 80a | 3931 MT Woudenberg | The Netherlands
cs@medicom-eu.com | +31 332 866 460