

H01038A
10.001D-10**Shoecovers without sole 100% PP**
Designed with recycled-origin materials
Personal Protective Equipment, Cat. I & Medical Device Class I

FEATURES & BENEFITS

- The hydrophobic spunbond of the shoe cover is made from **90% recycled-origin material**.
- Their anatomical shapes allow them to adapt to all types of shoes and foot shapes, particularly when used without shoes in medical environments.

Recommended uses

- In industrial environments on dry floors
- In medical environments



PRODUCT DETAILS

Product name :	Kolmi®
Product type:	Single-use, non-sterile
Composition:	Spunbond polypropylene Embossed polyethylene film
Spunbond weight :	30 g/m ²
Manufacturing :	Gluing
Size :	160 x 380 mm (Size L)
Color :	White
Unit weight:	4,8 g
Packaging :	Case of 500 pcs (Dispenser case or with protective bag)
Origin :	France

CERTIFICATIONS & STANDARDS

In compliance with the requirements of Regulation (EU) 2016/425 on Personal Protective Equipment.
Meet requirements of Regulation (EU) 2017/745 on Medical Devices.
Meet the requirements of the EN ISO 13688:2013+A1:2021 standard.
Manufacturing site: ISO 9001 and ISO 13485 certified.

MICROBIOLOGICAL INFORMATION

Evaluation of initial microbial contamination according to standards ISO 11737:2018 + A1:2021 upon request.
Additional microbiological tests: ASR, E. coli, staphylococci, available upon request.

H01038A
10.001D-10

Shoecovers without sole 100% PP
Designed with recycled-origin materials
Personal Protective Equipment, Cat. I & Medical Device Class I



MECHANICAL RESISTANCE DATAS (Cetelor Laboratory)

Tests on the polypropylene spunbond nonwoven	Standards	Results
Resistance to water penetration – increasing hydrostatic pressure	EN ISO 811:2018	13.4 cm H ₂ O
Abrasion resistance	EN ISO 12947-2:2017	6 000 cycles
Tensile strength	EN ISO 9073-3:2023	Machine direction: 70 N Cross direction: 55 N
Elongation at maximum force	EN ISO 9073-3:2023	Machine direction: 47 % Cross direction: 55 %
pH	EN ISO 3071:2020	6.0

NORMAL STORAGE CONDITIONS

Must not be exposed to moisture or sunlight, and should be stored in a cool or room temperature environment.

Product shelf life: 5 years.

LOGISTICS



**Outer case
spec.**

Item	Packaging Type	Qty pcs	Size (mm)	Gross weight (kg)	QTY/pallet
10.001D-10	Carton avec sache	500	400 x 400 x 310	2,7	36
H01038A	Carton distributeur	500	400 x 400 x 375	2,7	36



EAN Codes

Item	Outer case
10.001D-10	3 701 407 400 055
H01038A	3 701 407 400 048



Medicom SAS

Bld de la chanterie | 49124 Saint Barthélemy d'Anjou – France

T : 33 (0) 241 963 434 | F : + 33 (0) 241 963 453 | sales.support.pro.eu@medicom.com | www.medicom-eu.com